

**INTENDED USE**

The Zearalenone Strip Test Kit is a competitive chromatographic immunoassay for the semi-quantitative detection of the presence of Zearalenone in grains, feed materials, part of ready-product feed, spices. The assay could be used for on-site testing.

Detection Range: 5  $\mu$ g/kg - 50  $\mu$ g/kg (  $\mu$ g/kg = ppb)

**PRINCIPLE**

The Zearalenone Strip Test Kit is based on competitive lateral flow immunochromatographic assay. A test strip is placed inside a plastic cassette with a testing window. The specific antibodies for Zearalenone are labeled to the test system. The aflatoxins in specimen and the analytes complex pre-coated in the test strip would competitively bind to the antibodies. Result will be determined by if the colored red lines appears in the test zone.

**STORAGE AND STABILITY**

The kit should be stored at 2-8°C. The kit will be valid in 12 months. The test are stable through the expiration date printed on the label. **DO NOT FREEZE.** Do not use the kit beyond the expiration date.

**PRECAUTIONS**

- For best results, please strictly adhere to these instructions.
- Do not touch the membrane area of the test strips.
- Do not use the test beyond its expiration date marked on the foil pouch.
- The components in this kit have been quality control tested as standard batch unit. Do not mix components from different lot numbers.

**MATERIALS****Materials Provided**

- 8 dipsticks and 8 microwells sealed in one canister. 12 canisters each box.
- Assay buffer (2 bottle/box)
- Package insert for use

**Materials Required But Not Provided**

- Balance (sensitivity: 0.1g)
- Micropipettes (20-200 $\mu$ L, 100-1000 $\mu$ L)
- Centrifuge tube (2mL, 50mL)
- Centrifuge
- Timer

Homepage: [www.jgbiotech.com](http://www.jgbiotech.com)


**Solutions (according to sample extraction )**

- Sample extraction solution (50% ethanol): 50mL ethanol + 50mL deionized or distilled water

**SAMPLE PREPARATION AND EXTRACTION****Cut-offs: 5ppb, 10ppb, 20ppb, 30ppb, 50ppb**

- Grind a representative sample to the particle size of fine instant coffee and pass through a 20-mesh screen. Weigh out 5.0g of sample into a 50mL centrifuge tube.
- Add 12ml of Sample Extraction Solution (50% ethanol) into the tube and screw the cap. Shake for 3minutes.
- Let it stand or centrifuge for 5 minutes at 4000r/min.
- Transfer 30  $\mu$  L of supernatant layer of the extract into a new tube.
- Add assay buffer into the new tube for dilution according to the cut-off requirement listed in the following table. The diluted solution is the final Sample Extract for applying the assay.

Cut-off (ppb)	5	10	20	30	50
Assay Buffer( $\mu$ L)	150	400	850	1280	2150

**Cut-off: 2ppb**

- Grind a representative sample to the particle size of fine instant coffee and pass through a 20-mesh screen. Weigh out 5.0g of sample into a 50mL centrifuge tube.
- Add 12ml of Sample Extraction Solution (50% ethanol) into the tube and screw the cap. Shake for 3minutes.
- Let it stand or centrifuge for 5 minutes at 4000r/min.
- Transfer 60  $\mu$  L of supernatant layer of the extract into a new tube.
- Add 120  $\mu$  L of assay buffer into the new tube for dilution. The diluted solution is the final Sample Extract for applying the assay.

*NOTE: If the specimen is DDGS, corn germ cake, corn gluten cake, concentrate, corn fiber, corn bran, silage, and part ready-feed, please adjust the pH value of the supernatant*

Email: [info@jgbiotech.com](mailto:info@jgbiotech.com)

extract into pH 7.0-8.0 with 2.0mol/L of NaOH solution.

### **TEST PROCEDURE**

Read the product insert carefully before running the assay.

Allow the test kit, and specimen to room temperature (20-25°C) prior to testing.


1. Bring the kit to room temperature (20-25°C) before opening. Remove the test dipsticks and microwells from the canister and use it as soon as possible within 1 hour. Remain the rest strips and microwells in the canister and seal the cap immediately for future use.
2. Transfer 150 µL of the final Sample Extract into one microwell with a micropipette. Repeatedly suck and extrude the specimen until all red reagents are completely dissolved and start the timer.
3. Incubate the mixture for 3 minutes at room temperature (20-25°C). Then insert a dipstick into the microwell with the MAX end and start the timer.
4. Incubate for another 8 minutes. Take out the dipstick and interpret the result immediately according to the color intensities

### **INTERPRETATION OF RESULTS**

#### **NEGATIVE**

Both the test line and control line are visibly red.

#### **POSITIVE**



Colored C line is observed, no test line appears.

#### **INVALID**

Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test dipstick. If the problem persists,

discontinue using the test kit immediately and contact your local distributor.

### **QUALITY CONTROL**

Internal procedural controls are included in the test. A colored line appearing in the control region (C) is an internal procedural control. It confirms sufficient specimen volume and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

### **LIMITATIONS**

The Zearalenone Strip Test Kit is a useful tool offering a rapid and accurate screening test, especially for on-site mycotoxins testing. It provides a quick on-site detection method to detect Zearalenone in grains, feed, etc. If you want a quantitative result, it is suggested to apply other method such as ELISA or HPLC in practice.



J&G Biotech Ltd, (Reg. No.: 08419172)  
326 Cleveland Road, London,  
England E18 2AN, United Kingdom